

TABLE OF CONTENTS



OPERATING TABLES

ADDITIONAL ACCESSORIES FOR OPERATING TABLES

OBSTETRICS & GYNECOLOGY

PATIENT TRANSPORTERS

MULTIPURPOSE MEDICAL CHAIRS

MEDICAL FURNITURE

NOTATIONS







Manual control









Foot control







OM-DELTA-PLUS-02

GENERAL SURGERY OPERATING TABLE

- Electrohydraulic table with backup mechano-hydraulic drive
- Programmable table top positions










Used with the following KPP sets:

01 02 03 04 05 06 07 08 09 10 11 12 13 15 16 17 18 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43





Functions	OM-DELTA-PLUS-02
Table top height	
Table top lateral tilts	
Table top longitudinal tilts	
Table top zero position	
Table top longitudinal shift	
Dorsal section tilts	
Head section tilts	
Leg section tilts	®
Kidney bridge lift	
Dorsal section bending	<u> </u>
Extensible support shift	



Accessories







Foot control panel

4 I Operating tables I 5

OK-BETA PLUS UNIVERSAL OPERATING TABLE OUK-02-1

- Electrohydraulic table with backup mechano-hydraulic drive
- Programmable table top positions




Used with the following KPP sets:

01 02 03 04 05 06 07 08 09 10 11 12 13 14 15 16 17 18 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43





Functions	OK-BETA PLUS
Table top height	
Table top lateral tilts	
Table top longitudinal tilts	
Table top zero position	
Table top longitudinal shift	
Dorsal section tilts	
Head section tilts	®
Leg section tilts	<u> </u>
Dorsal section bending	
Extensible support shift	



Accessories







Pedal control unit

Operating tables | 7 6 | Operating tables

OUK-02-2 UNIVERSAL OPERATING TABLE

- Unique solution for thoracic surgery
- Electrohydraulic table with backup mechano-hydraulic drive
- The surgeon has the greatest possible access to the operating field










Used with the following KPP sets:

01 07 10 15 16 17 18 26 27 29 30 31 32 33 34 36 37 38 39 40 41 43





Functions	OUK-02-2
Table top height	
Table top lateral tilts	
Table top longitudinal tilts	
Table top zero positioning	
Table top longitudinal shift	
Dorsal section tilts	
Head section tilts	®
Leg section tilts	©
Kidney bridge lift	
Extensible support shift	



Accessories





Wireless remote control

Pedal control unit

8 | Operating tables | 9

OM-DELTA

GENERAL SURGERY OPERATING TABLES



- Electrohydraulically adjustable general surgery operating tables
- Programmable table top positions










Used with the following KPP sets:

01 02 03 04 05 06 07 08 09 10 11 12 13 15 16 17 18 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43





Functions	OM-DELTA-01	OM-DELTA-02
Table top height		
Table top lateral tilts		
Table top longitudinal tilts		
Table top zero position		
Table top longitudinal shift		
Dorsal section tilts		
Head section tilts	©	®
Leg section tilts	©	©
Kidney bridge lift		
Dorsal section bending	NO	
Extensible support shift		



Accessories



Operating tables I 11 control

10 | Operating tables





• New generation of high-functioning electromechanical tables for all types of surgical procedures

Kidney bridge

(OK-TETA PLUS-01)

• Programmable table top positions





sections without clamping

screws



Visco elastic memory

foam mattress



Dorsal section bending (OK-TETA PLUS)



Multipurpose rails made of stainless steel on all table top sections

> Pedal for activating the table to operating mode

Extensible supports to fix the table


> Base casing made of extra strong composite material





made of stainless steel






01 02 03 04 05 06 07 08 09 10 11 12 13 15 16 17 18 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43

Double castors Ø 100 mm





Functions	OK-TETA PLUS	OK-TETA PLUS-01
Table top height		<u> </u>
Table top lateral tilts		
Table top longitudinal tilts		
Table top zero position		
Table top longitudinal shift		
Dorsal section tilts		
Head section tilts	©	
Leg section tilts	<u> </u>	©
Kidney bridge lift	NO	
Dorsal section bending		NO
Extensible support shift	4	



Accessories





control

Pedal control panel

12 | Operating tables Operating tables | 13





 New generation of high-functioning electromechanical tables for all types of surgical procedures






01 02 03 04 05 06 07 08 09 10 11 12 13 14 15 16 17 18 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43

Double castors Ø 100 mm





Functions	OK-TETA	OK-TETA-01
Table top height		
Table top lateral tilts		
Table top longitudinal tilts		
Table top zero position		
Table top longitudinal shift		
Dorsal section tilts		
Head section tilts	©	©
Leg section tilts	©	©
Dorsal section bending	NO	
Extensible support shift	4	42



14 | Operating tables | 15



I SARA MORE

• New generation of high-functioning electromechanical tables for all types of surgical procedures




Used with the following KPP sets:

01 02 03 04 05 06 07 08 09 10 11 12 13 14 15 16 17 18 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43





Functions	OK-BETA	OK-BETA-01	OK-BETA-02	OK-BETA-03
Table top height				
Lateral tilts				
Table top longitudinal tilts				
Table top zero position				
Table top longitudinal shift	NO	NO		
Dorsal section tilts	•	•	•	•
Head section tilts	•	•	•	•
Leg section tilts	•	•	•	•
Dorsal section bending	NO		NO	
Extensible support shift	4	<i>\(\pi\)</i>	\(\pi\)	4




Operating tables I 17




MEDIN ALFA OPERATING TABLE OU-01K




Functions	OU-01K	OU-01K-01	OU-01K-02	OU-01K-03
Table top height	4	Ø	4	4
Lateral tilts	Ø	100	4	4
Table top longitudinal tilts	B	100	1	4
Table top longitudinal shift	NO	NO		
Dorsal section tilts	•	•	®	(
Head section tilts	•	•	®	(
Leg section tilts	•	•	•	(h
Dorsal section bending	NO		NO	
Extensible support shift	12	<i>\(\opi\)</i>		4



- Universal operating tables with mechano-hydraulic adjustment for all types of surgical procedures
- Ideal for the most challenging work environments



Molded polyurethane foam cushions. Made of latex-free, fire-resistant, antistatic material





Pedal for activating the table to operating mode

Base casing made of extra strong composite material

Extensible supports to fix the table



фмедин"



Double castors

Ø 100 mm


made of stainless steel

Height/tilt adjustment

pedal







01 02 03 04 05 06 07 08 09 10 11 12 13 14 15 16 17 18 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43

Table top tilt control

18 | Operating tables Operating tables I 19





Mechano-hydraulically adjustable universal operating table










Used with the following KPP sets:

01 05 07 08 10 12 14 15 16 17 18 26 27 29 30 31 32 33 34 35 36 37 38 39 40 41 43





Functions	OUK-01
Table top height	
Table top lateral tilts	©
Table top longitudinal tilts	©
Dorsal section tilts	©
Head section tilts	©
Leg section tilts	©
Extensible support shift	



20 I Operating tables 1 21



OPERATING TABLES

antistatic material



- Mechano-hydraulically adjustable general surgery operating tables
- Easy to move due to antistatic double castors with a central-locking brake






Table top tilt control



Column protective covers and base casing made of extra strong composite material




Height adjustment pedal



Central brake and

Double castors Ø 150 mm

motion direction fixation pedal




01 05 07 08 10 12 14 15 16 17 18 26 27 29 30 31 32 33 34 35 36 37 38 39 40 41 43





Functions	OK-GAMMA MOBILE 01	OK-GAMMA MOBILE 02	OK-GAMMA MOBILE 03
Table top height	4	100	1
Table top lateral tilts			NO
Table top longitudinal tilts	•	(©
Dorsal section tilts	®	•	(
Head section tilts	®	(©
Leg section tilts	•	•	©
Dorsal section bending		NO	NO



22 | Operating tables Operating tables I 23





- Mechano-hydraulically adjustable general surgery operating tables
- Optimal solution for patient preparation and recovery, minor surgical procedures, day surgery










Used with the following KPP sets:

01 05 07 08 10 12 14 15 16 17 18 26 27 29 30 31 32 33 34 35 36 37 38 39 40 41 43





Functions	OK-GAMMA 01	OK-GAMMA 02	OK-GAMMA 03
Table top height	B	4	1
Table top lateral tilts		NO	
Table top longitudinal tilts	•	•	(
Dorsal section tilts	®	®	©
Head section tilts	(b)	((
Leg section tilts	(•	@
Dorsal section bending		NO	NO




Operating tables I 25







Functions	MEDIN-SAFIS	MEDIN-SAFIS 02
Table top height		
Table top lateral tilts		
Table top longitudinal tilts		
Table top longitudinal shift		NO
Table top lateral shift		NO
Table top zero position		
Dorsal section tilts	NO	NO
Switching the table to the fixed position	4	



Accessories



for accessories mounting




Cross bar with head rest



Table side bar extension



Arm protection





Injection plate



Patient restraint strap



26 | Operating tables Operating tables I 27

MEDIN-SAFIS 03

RADIOGRAPHIC SURGERY TABLE

- Electromechanically adjustable radiographic operating table
- Two-speed longitudinal-lateral shift mode
- Programmable table top positions




15 16 18 26 30 31 32 33 37 38 39 40 41 43






Functions	MEDIN-SAFIS 03
Table top height	
Table top lateral tilts	
Table top longitudinal tilts	
Table top longitudinal shift	
Table top lateral shift	
Table top zero position	
Dorsal section tilts	©
Switching the table to the fixed position	♥ •



















Accessories



Leg section



Surgical leg holders (recommended)



side holder



Instrument table



for infusions



Injection plate



Patient restraint strap



Lateral support

28 | Operating tables Operating tables I 29





- Portable folding operating table
- Simple design allows the general surgery portable folding operating table to be assembled in just a few movements, so it is immediately ready for use in the field and extreme environments

• Transported in two special containers



Head section with mechanical fixation

Molded polyurethane foam cushions. Made of latex-free, fire-resistant, antistatic material



Removable handles for table top tilt adjustment

> Base frame mad of stainless steel

Screw supports to set the table into the working position



Multipurpose rails made of stainless steel

Table top parts and accessories are stored in the lodgements of the first container (dimensions: 1,200x800x510 mm)




Table lifting/ lowering pedal



Base frame is fixed on the 2nd container





01 05 07 10 12 15 16 17 18 26 27 29 30 31 32 33 34 35 36 37 38 39 40 41 43


Double castors Ø 100 mm


(equipped with brakes)






Functions	OK-OMEGA 01	OK-OMEGA 02
Table top height		42
Table top lateral tilts	(®
Table top longitudinal tilts	(o
Dorsal section tilts	(•
Head section tilts	(®
Leg section tilts	(o
Dorsal section bending	NO	











30 | Operating tables Operating tables | 31



OM-SIGMA

OPERATING TABLE SYSTEMS WITH INTERCHANGEABLE TABLE TOPS



- The unique all-in-one system that meets all up-to-date requirements and offers infinite possibilities, thus making the entire process in the surgical suite optimal
- It is a combination system that includes an operating table and a patient transporter
- Preserves a sterile environment and increases the patient capacity of the surgical suite
- No step of transferring the patient to the patient transporters, thereby reducing the risk of post-operative injury
- Three types of table tops allow all types of surgical procedures to be performed on the table
- Programmable table top positions












Table top transport trolleys



















Visco elastic memory Head section with dual tilt function

Column protective covers and base casing made of stainless steel



The flat base of the column provides more space for operating the C-arm X-ray machine in the lower position of the table

Emergency stop button

32 | Operating tables | 33



OPERATING TABLE SYSTEMS WITH INTERCHANGEABLE TABLE TOPS


Interchangeable table tops

- The OM-SIGMA table can be used with three types of table tops: universal, special, and diagnostic
- Three types of table tops allow the operating table systems to be used for different surgical specialties

Universal table top



It is intended for all types of surgical procedures. The table top is X-ray-transparent. The table top has side rails for attaching accessories.




Special table top



It is intended for performing traumatological and orthopedic surgical interventions on the lower limbs and shoulder joint.

The surgeon has a comfortable working position and good access to the operating field.


Table top sections can be detached while the patient is on the table.



Diagnostic table top



It is made of carbon fiber. It is designed for a wide range of interventional diagnostic and therapeutic radiosurgical procedures using the C-arm X-ray machine, as well as part of the angiographic unit.



OM-SIGMA

OPERATING TABLE SYSTEMS WITH INTERCHANGEABLE TABLE TOPS



Operating table column

- The guide column includes an electrohydraulic drive that controls all table top movements: lifting-lowering, lateral left-right tilting, Trendelenburg and reverse Trendelenburg positioning
- The batteries are inside the column
- The main positions of the table can be controlled remotely using a wireless remote control or a backup control panel on the column







Column (moved by the transporter)

Mobile column

Stationary column

Column type	Table top type	Operating table system with interchangeable table tops	Safe working load	Maximum table load
	universal	OM-SIGMA-06	250 kg	450 kg
Column (moved by the transporter)	special	OM-SIGMA-07	250 kg	450 kg
	diagnostic	OM-SIGMA-08	190 kg	250 kg
	universal	OM-SIGMA-09	250 kg	450 kg
Mobile	special	OM-SIGMA-10	250 kg	450 kg
	diagnostic	OM-SIGMA-11	190 kg	250 kg
	universal	OM-SIGMA-12	250 kg	450 kg
Stationary	special	OM-SIGMA-13	250 kg	450 kg
	diagnostic	OM-SIGMA-14	190 kg	250 kg

Used with the following KPP sets:

for the universal table top

01 02 03 04 05 06 07 08 09 10 11 12 13 15 16 17 18 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43

for the special table top


01 02 03 04 05 06 07 08 15 16 18 26 27 29 30 31 32 34 36 37 38 39 40 41 43


for the diagnostic table top

01 15 16 18 26 30 31 32 37 38 39 40 41 43

Operating tables I 35

	SAFIS	SAFIS-03	OK-DELTA PLUS-02	OK-DELTA	OK-BETA PLUS	OUK-02-2	OK-TETA PLUS	OK-BETA	OK-TETA	MEDIN ALFA	OK- EPSILON	OK- OMEGA	OK- GAMMA	OK-GAMMA MOBILE
KPP-01			+	+	+	+	+	+	+	+	+	+	+	+
KPP-02			+	+	+		+	+	+	+				
KPP-03			+	+	+		+	+	+	+				
KPP-04			+	+	+		+	+	+	+				
KPP-05			+	+	+		+	+	+	+	+	+	+	+
KPP-06			+	+	+		+	+	+	+				
KPP-07			+	+	+	+	+	+	+	+	+	+	+	+
KPP-08			+	+	+		+	+	+	+	+		+	+
KPP-09			+	+	+		+	+	+	+				
KPP-10			+	+	+	+	+	+	+	+	+	+	+	+
KPP-11			+	+	+		+	+	+	+				
KPP-12			+	+	+		+	+	+	+	+	+	+	+
KPP-13			+	+	+		+	+	+	+				
KPP-14					+			+	+	+	+		+	+
KPP-15	+	+	+	+	+	+	+	+	+	+	+	+	+	+
KPP-16	+	+	+	+	+	+	+	+	+	+	+	+	+	+
KPP-17			+	+	+	+	+	+	+	+	+	+	+	+
KPP-18	+	+	+	+	+	+	+	+	+	+	+	+	+	+
KPP-26	+	+	+	+	+	+	+	+	+	+	+	+	+	+
KPP-27			+	+	+	+	+	+	+	+	+	+	+	+
KPP-28			+	+	+		+	+	+	+				
KPP-29			+	+	+	+	+	+	+	+	+	+	+	+
KPP-30	+	+	+	+	+	+	+	+	+	+	+	+	+	+
KPP-31	+	+	+	+	+	+	+	+	+	+	+	+	+	+
KPP-32	+	+	+	+	+	+	+	+	+	+	+	+	+	+
KPP-33		+	+	+	+	+	+	+	+	+	+	+	+	+
KPP-34			+	+	+	+	+	+	+	+	+	+	+	+
KPP-35			+	+	+		+	+	+	+	+	+	+	+
KPP-36			+	+	+	+	+	+	+	+	+	+	+	+
KPP-37	+	+	+	+	+	+	+	+	+	+	+	+	+	+
KPP-38	+	+	+	+	+	+	+	+	+	+	+	+	+	+
KPP-39	+	+	+	+	+	+	+	+	+	+	+	+	+	+
KPP-40	+	+	+	+	+	+	+	+	+	+	+	+	+	+
KPP-41	+	+	+	+	+	+	+	+	+	+	+	+	+	+
KPP-42			+	+	+		+	+	+	+				
KPP-43	+	+	+	+	+	+	+	+	+	+	+	+	+	+
















1. Anesthetic screen	1 pc.
2. Medical stand for infusions (SHDV	-5) 1 pc.
3. Injection plate	2 pcs.
4. Patient restraint strap	1 pc.
5. Hand strap	2 pcs.
6. X-ray cassette holder	1 pc.
7. Lateral support	2 pcs





















1 pc.

1 pc.

1 pc.

KPP-02 set includes the following:

1. Orthopedic section

10. Trolley





38 | Additional accessories for operating tables Additional accessories for operating tables 1 39





SET OF ADDITIONAL ACCESSORIES FOR ORTHOPEDICS AND TRAUMATOLOGY SURGERIES ON LOWER LEG AND KNEE






Addition to the basic set.
Used with KPP-02 set only











1. Lifting assembly with perineal support (thigh support)	1 pc.
2. Extension for orthopedic section	1 pc.
3. Leg support with radial clamp	1 pc.
4. Rotary lever	1 pc.
5. Tension bracket holder	1 pc.





KPP-04 set includes the following:

Universal swivel popliteal support with hinged stabilization of joint

1 pc.

1 pc.

2. Tension bräcket holder






40 I Additional accessories for operating tables I 41





SET OF ADDITIONAL ACCESSORIES FOR ORTHOPEDICS AND TRAUMATOLOGY SURGERIES (FOR PLASTER BANDAGING IN PELVIC ZONE AND LOWER LIMBS)





















5. Hand stabilizing cuff







KPP-06 set includes the following:

1. Assembly for plaster bandaging in pelvic zone 1 pc. 2. Cross bar stand for leg suspending 1 pc.







42 | Additional accessories for operating tables Additional accessories for operating tables 1 43

1 pc.







2 pcs.











KPP-07 set includes the following:

١.	Popliteal support with adjustable thigh holder	I pc.
2.	Gear unit with swivel clamp of popliteal support	1 pc.
3.	Operation crank	1 pc.
1	Crank holder	1 no





KPP-08

SET OF ADDITIONAL ACCESSORIES FOR ORTHOPEDICS AND TRAUMATOLOGY SURGERIES (FOR BONE TRACTION)







KPP-08 set includes the following:

1. Rotation unit 1 pc. 2. Tensioner 1 pc.

Used as an addition to KPP-02, KPP-03, KPP-04, KPP-05 sets






















1.	Table mounting adapter	1 pc.
2.	Hinge-type holder	1 pc.
3.	Dorsal support with cushion	1 pc.
4.	Connecting element	1 pc.
	Ball-type connecting element	1 pc.
5.	Cranial clamp	1 pc.
6.	Pins to cranial clamp (4 types)	3 pcs each
7.	90 mm rest with support base	2 pcs.
8.	Adjustable head rest (left and right)	1 pair
	Set of head rest cushions	1 pair
9.	Ring-shaped head cushion assembly, 140 mm	1 pc.
10.	90 mm ring-shaped cushion with head support base	1 pc.
11.	Diametrical attachment	1 pc.





12. Radial serrated clamp

44 I Additional accessories for operating tables Additional accessories for operating tables 1 45







KPP-10 set includes the following:

1. Leg support with radial serrated clamp

2. Bowl (6.5 L) with side holder

2 pcs. 1 pc.







KPP-11 set includes the following:

Leg support with radial serrated clamp
 Bowl (6.5 L) with side holder and drainage system
 pc.

3. Extension section 1 pc.





Additional accessories for operating tables I 47



SET OF ADDITIONAL ACCESSORIES FOR ENT AND OPHTHALMOLOGY SURGERIES

























KPP-12 set includes the following:

Support bracket		1 pc.
Supporting leg with for head rests	n a mounting bar	1 pc.
3. Cross rail		1 pc.
4. Surgeon wrist rest		1 pc.
5. Stand		1 pc.
6. Ring-shaped cushi	on base	1 pc.
7. Padding disk		1 pc.
8. Ring-shaped cushi	on	1 pc.
9. Table mounting ad	lapter	1 pc.
10. Surgeon arm rest		1 pc.
11. Adjustable U-shape (left and right head		1 pair
12. Set of cushions for (left and right head		1 pair











KPP-13 set includes the following:

1. Roller	1 pc.
2. Stand 1	1 pc.
3. Stand 2	1 pc.
4. Knee section (with in-built clamps)	1 set
5. Radial serrated clamp	2 pcs.






SET OF ADDITIONAL ACCESSORIES FOR SURGERIES REQUIRING BODY EXTREMITIES POSITIONING (ELEVATION AND STABILIZING OF LUMBAR, PELVIC OR CHEST ZONES)









KPP-14 set includes the following:

1.	Lifting assembly	1 pc.
2.	Pad of lifting zone	1 pc.
3.	Central section pad	1 pc.
4.	Plate	1 pc.
5.	Crank handle	1 pc.



KPP-16

SET OF ADDITIONAL ACCESSORIES FOR SURGERIES REQUIRING BODY EXTREMITIES POSITIONING (LOWER POSITIONING OF ARM)





KPP-17

KPP-16 set includes the following:

1. Radial serrated clamp	1 pc
2. Arm rest	1 pc
3. Support bracket	1 pc




KPP-15

SET OF ADDITIONAL ACCESSORIES FOR SURGERIES REQUIRING BODY EXTREMITIES POSITIONING (UPPER POSITIONING OF ARM)



SET OF ADDITIONAL ACCESSORIES FOR SURGERIES REQUIRING BODY EXTREMITIES POSITIONING (LATERAL POSITIONING OF ARM)





KPP-15 set includes the following:

1. Radial serrated clamp 1 pc. 2. Arm rest with supporting bar 1 pc.



KPP-17 set includes the following:

1.	Supporting leg with sliding block	2 pcs.
2.	Frame with ball joint	1 pc.
3.	Pad	1 pc.


50 | Additional accessories for operating tables Additional accessories for operating tables | 51



SET OF ADDITIONAL ACCESSORIES FOR SURGERIES REQUIRING BODY EXTREMITIES POSITIONING (LEG POSITIONING ACCORDING TO GOEPEL)








Leg support	1 pc.
Radial serrated clamn	1 nc



KPP-27

SET OF ADDITIONAL ACCESSORIES FOR SURGERIES REQUIRING OPERATING TABLE EXTENSION (FOR TABLE TOP WIDENING)





KPP-27 set includes the following:

1. Support 1 pc. 2. Pad 1 pc.



KPP-26

SET OF ADDITIONAL ACCESSORIES FOR SURGERIES REQUIRING HEAD POSITIONING (RING-SHAPED HEAD CUSHIONS)



KPP-28


SET OF ADDITIONAL ACCESSORIES FOR SURGERIES REQUIRING OPERATING TABLE EXTENSION (FOR TABLE TOP EXTENDING)





KPP-26 set includes the following:

1.	Ring-shaped head cushion, Ø 205 mm	1 pc.
2.	Ring-shaped head cushion, Ø 140 mm	1 pc.
3	Ring-shaped head cushion Ø 90 mm	1 nc



KPP-28 set includes the following:

1. Frame 1 pc. 2. Pad 1 pc. 3. Clamping screw 2 pcs.

52 | Additional accessories for operating tables Additional accessories for operating tables 1 53



SET OF ADDITIONAL ACCESSORIES FOR SURGERIES REQUIRING OPERATING TABLE EXTENSION (FOR TABLE SIDE BAR ELONGATION)







1.	160 mm extension bar	1 pc.
2.	260 mm extension bar	1 pc.
3	460 mm extension har	1 nc



KPP-31

SET OF ADDITIONAL ACCESSORIES FOR SURGERIES REQUIRING OPERATING TABLE EXTENSION (STEP PLATES TO OPERATING TABLE)





KPP-31 set includes the following:

	300 mm step plate	1 pc
)	220 mm step plate	1 pc
2	130 mm sten plate	1 nc



KPP-32

SET FOR ACCESSORIES STORAGE (MOBILE STAND)





KPP-30

SET OF ADDITIONAL ACCESSORIES FOR SURGERIES REQUIRING OPERATING TABLE EXTENSION (FOR EXTENSION OF ANESTHETIC SCREEN)





KPP-30 set includes the following:

1. Extension pieces to anesthetic screen 2 pcs.





Used with the following KPP sets:

01 02 03 04 05 06 07 08 10 11 12 13 15 16 17 18

KPP-32 set includes the following:

1. Stand 1 pc. 2. Bottom tray



54 | Additional accessories for operating tables



SET OF SURGICAL LEG HOLDERS






KPP-34 LIFTING BRIDGE FOR SPINAL SURGERY











KPP-33 set includes the following:

1. Left surgical leg holder 1 pc. 2. Right surgical leg holder 1 pc.

















KPP-34 set includes the following:

1. Bridge	1 pc.
2. Pad	2 pcs.
3. Head pad	1 pc.
4. Arm pad	2 pcs.
5. Knee pad	1 pc.
6. Hemispheric cushion	1 pc.
7. Crank handle	1 pc.




56 | Additional accessories for operating tables Additional accessories for operating tables 1 57





KPP-35 set includes the following:

1. Dorsal section	1 pc.
2. Right shoulder section	1 pc.
3. Left shoulder section	1 pc.
4. Pad	1 pc.
5. Lateral support	1 pc.
6. Leg support with radial serrated clamp	1 pc.



KPP-36 set includes the following:

	Base plate	1 pc.
)	Guide rail	1 pc.
3.	Sliding member	1 pc.
1.	Splint	1 pc.
).	Clamp	2 pcs
).	Bracket	1 pc.
7.	Support	1 pc.
3.	Leg support	1 pc.

58 I Additional accessories for operating tables I 59







Handle for easy transferring













2. Patient lateral positioning cushion KPP-38



3. Spinal surgery cushion KPP-39



4. Patient leg cushion KPP-40



KPP-41 set includes the following:

Rigid frame
 Movable fabric cover

1 pc. 1 pc.

Dimensions (L/W/H): 1,495x490x7 mm






60 | Additional accessories for operating tables | 61








KPP-42 set includes the following:

1. Carbon fiber panel

2. Mattress







KPP-43 set includes the following:

1. Quick-release clamp


2. Arm rest





- EARN MORE
- Provides high level of comfort, safety and ergonomics for the parturient woman and midwifery team
- Comfortable patient transportation within the healthcare facility
- Can be used in bed and chair positions with extensible leg section

Removable anti-bedsore mattress and covers made of waterproof and breathable fabric





Feet rests



Remote controls









Trolley for accessories







Functions	KA-01	KA-02	KA-03	KA-04
Table top height				4
Longitudinal tilt		<u>M</u>	•	(
Dorsal section tilt				•
Leg section height		NO	NO	NO
Leg section tilt		NO	NO	NO
















Accessories



for infusions











for magazines Sci

64 | Obstetrics & Gynecology | 65

SR-01 OBSTETRIC AID TABLE

EARN MORE

- Provides high level of comfort for the parturient woman and midwifery team
- Comfortable patient transportation within the healthcare facility
- The table can be used in chair or bed position













Functions	SR-01
Table top height	
Seat tilts	
Backrest tilts	
Head section tilts	©
Height of the table for the newborn	©





Chair position

Accessories



for infusions

Injustion plate





Injection plate Bowl with drainage

66 Obstetrics & Gynecology 67

KGM-4P GYNECOLOGICAL CHAIR

- Chair with electrical adjustment of height, backrest tilts, seat tilts, arm rests with foot rests
- Various leg support systems can be used: arm rests with foot rests, arm rests with Goepel leg holders
- Initial low seating position and the arm rest with foot rest that can be tilted backwards to facilitate patient positioning make the chair comfortable and safe to use for people with disabilities
- It can be used in chair or bed position. Enhanced functionality up to operating table level










Upholstery colors








Functions	KGM-4P
Seat height	
Seat tilts	
Backrest tilts	
Height of arm rests with foot rests	
















Table position

Accessories







Injection plate



drainage



Leg section



Lightweight leg Step plate



Side rails (left/right)



Arm rest with Goepel leg holder (left/right)



Cover for Goepel leg holder



Colposcope holder





control unit



remote control

section



Obstetrics & Gynecology | 69





• Designed for examinations, manipulations, minor and major gynecological surgeries, urological procedures, colposcopy




3096 3022 1549 1017 1089 1210 1044 6156 6279 6001 6099 5154 5141 5118 3036 7000 9011 9001







Functions	KGM-3P	KGM-2P	KGM-1
Seat height	<i>€</i> 2	€2	NO
Seat tilts	<i>€</i> 2	€2	\$
Backrest tilts	\S	NO	NO
















Table position (KGM-3P)

Accessories







leg holder









Colposcope holder

drainage





Lightweight

leg section*

* Only for KGM-3P Obstetrics & Gynecology | 71

KGM-2, 4 GYNECOLOGICAL CHAIRS



- Designed for examinations, manipulations, minor and major gynecological surgeries, urological procedures, colposcopy
- It can be used in chair or bed position. Enhanced functionality up to operating table level




Upholstery colors

3096	3022	1549	1017	1089	1210	1044	6156	6279	6001	6099	5154	5141	5118	3036	7000	9011	9001







Functions	KGM-2	KGM-4
Seat tilt	<i>€</i> 2	(
Backrest tilt	<i>₹</i>	(















Table position

Accessories











leg section



drainage

Colposcope holder

Side rails (left/right) Cover for

Goepel leg holder





Pedal control unit* Disposable sheet holder

* Only for KGM-2 Obstetrics & Gynecology | 73





- Comfortable ergonomic design and dimensions
- Smooth and noiseless bed tilt adjustment by means of gas spring
- The design allows the newborn bed to be placed directly at the mother's bedside




 $[\]ensuremath{^{\star}}$ The image of the mattress is just promotional

	KGM-4P	KGM-3P	KGM-2P	KGM-1	KGM-2	KGM-4
Electric height adjustment	+	+	+			
Electric seat tilt adjustment	+	+	+	+	+	
Electric backrest tilt adjustment	+	+			+	
Electric adjustment of arm rests with foot rests	+					
Seat tilt adjustment by means of gas springs						+
Backrest tilt adjustment by means of gas springs						+
Head rest adjustment	+	+	+	+	+	+
Horizontal positioning of backrest and seat	+	+			+	+
Leg section attachment is possible	+	+			+	+
Colposcope attachment is possible	+	+	+	+	+	+
Built-in socket for additional devices	+	+	+	+		
Paper towel holder	+	+	+	+	Р	Р
Rails (on the backrest) for accessories mounting	+	+	+	+	+	+
Rails (on the seat) for accessories mounting	+	+	+	+	+	+
Leg support height adjustment	+	+	+	+	+	+
Leg support tilt angle adjustment	+	+	+	+	+	+
Leg supports can be rotated	+	+	+	+	+	+
Withdrawable gynecological bowl (9 L)	+	+	+	+		
Withdrawable gynecological bowl (3.5 L)					+	+
Arm rests with foot rest	+					
Arm rests with Goepel leg holders	Р	+	+	+		
Goepel leg holder					+	+

P - Accessory + - Availability

74 Obstetrics & Gynecological chairs comparison 75









Medical stand for infusions

- Comfortable and maneuverable patient transporter with foot-operated hydraulic height and longitudinal tilt adjustment
- May be operated from both sides
- X-ray transparent table top with X-ray cassette tunnel
- Removable mattress dimensions: 1,950x640x80 mm























Central brake

and motion direction fixation control

Table top lifting/lowering

pedals (on both sides)

Accessories







Head rest

Patient rollboard holder



X-ray cassette

holder

Patient restraint



Injection plate

Ø 200 mm rubbered castors

(equipped with a central-locking brake)





Mesh basket

MEDIN-ELPIS 02 RESUSCITATION STRETCHER




Medical stand for infusions

- Comfortable and maneuverable patient transporter with foot-operated hydraulic height adjustment
- May be operated from both sides

• X-ray transparent table top with X-ray cassette tunnel













Accessories



Motion direction

fixation





holder



holder



strap





Multipurpose plate

Patient rollboard

TBP-01 PATIENT TRANSPORT STRETHERS

- Patient transporter with foot-operated hydraulic height adjustment
- Removable mattress dimensions: 1,870x600x80 mm



Accessories



Motion direction fixation



Head rest



Oxygen bottle holder



Chest strap/ leg strap




Patient rollboard

TBP-02 PATIENT TRANSPORT STRETHERS



- Patient transporter with foot-operated hydraulic height adjustment
- May be operated from both sides
- Mattress dimensions: 1,870x600x60 mm



Accessories



Motion direction fixation



Chest strap/ leg strap



for infusions

80 | Patient transporters Patient transporters | 81

TBPU PATIENT WASHING STRETCHERS



- Dedicated maneuverable patient transporter with foot-operated hydraulic height adjustment
- Ideal for caring for critically ill patients and patients with limited mobility



Accessories



Motion direction fixation

TB-01 PATIENT TRANSPORT STRETCHERS



- Lightweight fixed-height patient transporter with removable stretcher
- Designed for transfer and transportation of patients within healthcare facilities, as well as delivery to ambulance vehicles
- Mattress dimensions: 1,915x480x30 mm



Accessories




Chest strap/ leg strap


82 | Patient transporters | 83

TAP-01 PATIENT TRANSPORT STRETCHERS



- Maneuverable patient transporter with foot-operated hydraulic height adjustment
- Used in prosectoriums and morgues
- May be operated from both sides





KMMP

MULTIPURPOSE MOBILE MEDICAL CHAIRS



- Used for rehabilitation units, nursing centers, treatment rooms, cosmetology, ultrasound, and other rooms in medical facilities
- Highly stable, functional and maneuverable








Upholstery colors










Functions	KMMP-01	KMMP-02	KMMP-03
Chair height	4	4	
Additional simultaneous seat and backrest tilt (Help function)		NO	NO
Backrest longitudinal tilt	•	•	•
Leg section longitudinal tilt	•	®	<u> </u>
Leg section adjustment by patient			NO













Simultaneous backrest and leg section adjustment by means of gas springs



Convenient arm rest positions



Independent leg section adjustment by patient (KMMP-01, KMMP-02)



Height is adjusted by hydraulic drive (by pedal)

Accessories



Lateral body

supports



Head rest with lateral restraints





for dining



for infusions



Withdrawable foot plate (KMMP-02, KMMP-03)

KMM

MULTIPURPOSE MEDICAL CHAIRS



- Designed for blood donation centers, blood transfusion centers, treatment rooms, cosmetology, ultrasound, and other rooms in medical facilities
- Selection and setting the seat tilt angle in one of three locked positions (0, 7, 17°)









MEDICAL CABINETS

MEDICAL MODULE SYSTEMS



- Combination of medical cabinets integrated into a unified system to create optimal working environment for medical staff
- Designed for storage of medicines, pharmaceuticals, instruments, dressings, medical devices, and tools in medical centers, hospital departments, polyclinics, pharmacies



MEDICAL TALL CABINET UNITS

MEDICAL MODULE SYSTEMS



Medical tall cabinet unit for storing medicines (with shelves)

Dimensions (W/D/H):

- 300x600x2100 mm (707)
- 450x600x2100 mm (705)
- 600x600x2100 mm (701)



Medical tall cabinet unit for storing medicines (with shelves)

Dimensions (W/D/H):

- 900x600x2100 mm (708)



Medical tall cabinet unit for storing medicines for installing the in-built refrigerator (712) for installing built-in fridge and safe (713)

Dimensions (W/D/H): 600x600x2100mm



Medical tall cabinet unit for storing medicines (with shelves and mezzanine cabinet)

Dimensions (W/D/H): 900x600x2700mm (717)



In-built refrigerator (not included in the delivery set)



Withdrawable plate and in-built safe (dimensions of the safe (W/D/H): 550x480x270 mm)

2003 1015 6027 5012 4003 7035 9016



MEDICAL TALL CABINET UNITS

MEDICAL MODULE SYSTEMS



Medical tall cabinet unit for storing medicines (with withdrawable sections)

Dimensions (W/D/H): - 600x600x2100 mm (710)



Medical tall cabinet unit for storing medicines

Dimensions (W/D/H):

- 300x600x2100 mm (715)
- 600x600x2100 mm (711)



Catheter holders



Height-adjustable removable trays



Medical tall cabinet unit for storing catheters

Dimensions (W/D/H): - 600x600x2100 mm (741)





Medical tall cabinet unit for storing endoscopes (740) for storing endoscopes (glass door) (743)

Dimensions (W/D/H): - 600x600x2100 mm



Recirculator-irradiator



Simultaneous placement of 8 endoscopes



MEDICAL LOWER CABINET UNITS

MEDICAL MODULE SYSTEMS



Medical lower cabinet unit for storing tools and dressings (with shelves)

Dimensions (W/D/H):

- 450x600x850 mm (724)
- 600x600x850 mm (723)



Medical lower cabinet unit for storing tools and dressings (with shelves)

Dimensions (W/D/H):

- 900x600x850 mm (725)



Medical lower cabinet unit for medicines (with drawers)

Dimensions (W/D/H):

- 450x600x850 mm (728)
- 600x600x850 mm (727)



Medical lower cabinet unit for medicines (with drawers)

Dimensions (W/D/H):

- 900x600x850 mm (726)



Medical lower cabinet unit for wash basin installation

Dimensions (W/D/H):

- 450x600x850 mm (730)
- 430x370x180 mm



Medical lower cabinet unit for wash basin installation

Dimensions (W/D/H):

- 1200x600x850 mm (731)
- 900x600x850 mm (738)

Dimensions (W/D/H):

- 430x370x180 mm



Medical lower cabinet unit with slop hopper

Dimensions (W/D/H): - 600x600x850 mm (732)



Medical lower cabinet unit for installation of a gypsum sump

Dimensions (W/D/H):

- 1200x600x850 mm (734**)**

92 | Medical furniture | 93



MEDICAL LOWER CABINET UNITS

MEDICAL MODULE SYSTEMS



Medical lower cabinet unit for medical waste (with withdrawable sections)

Dimensions (W/D/H): - 300x600x850 mm (729)



Medical lower cabinet unit for medical waste (with withdrawable sections)

Dimensions (W/D/H):

- 600x600x850 mm (735)



Medical lower cabinet unit for medical waste (with two withdrawable sections)

TABLE TOPS. ACCESSORIES

MEDICAL WALL CABINET UNITS

Dimensions (W/D/H): - 600x600x850 mm (736)

MEDICAL MODULE SYSTEMS

MEDICAL MODULE SYSTEMS

Medical wall cabinet unit

Dimensions (W/D/H):

- 450x340x600 mm (720)

- 600x340x600 mm (721)

for dressings

(with shelves)

• Not included in the basic package, to be ordered additionally

• Table tops can be made with one or more (up to four)

can be integrated into one seamless element up to 2,400 mm long

• Table tops for side-by-side cabinets

• Made of stainless steel

set-in wash basins (430x370x180 mm)



Medical lower cabinet unit for storing medicines (for installing the in-built refrigerator) (the fridge is not included in the delivery set)

Dimensions (W/D/H): - 600x600x850 mm (737)

Medical wall cabinet

Dimensions (W/D/H):

900x340x600 mm (722)

unit for dressings

(with shelves)












SINGLE-DOOR, DIMENSIONS (W/D/H): 455X460X1,850 MM



MSH-1-01 (without safe) MSH-1-02 (with safe)



MSH-1-03



MSH-1-04

MSH-1-06



MSH-1-07

MSH-1-05

TWO-DOOR, DIMENSIONS (W/D/H): 900X460X1,850 MM



MSH-2-01 (without safe) MSH-2-02 (with safe)



MSH-2-03



MSH-2-04



MSH-2-05



MSH-2-06



MSH-2-07



94 | Medical furniture Medical furniture | 95





- Designed to equip the central sterilization department (CSD), hospital depots, and stockrooms of various departments
- Made of premium stainless steel
- 4 removable, height-adjustable shelves
- Simple and easy to assemble, shelves can be rearranged without tools



Dimensions (versions): length: 830/1030/1230 mm width: 600/700 mm height: 1750 mm Ø 125 mm castors



rigid construction

2

Heavy load she



Pins for shelves and rails



Castors with brakes



SMS/SMSP

MEDICAL STATIONARY/SIDE RACKS

- · Designed to equip the central sterilization department (CSD), hospital depots, and stockrooms of various departments
- Made of premium stainless steel
- 6 removable, height-adjustable shelves
- Simple and easy to assemble, shelves can be rearranged without tools



Dimensions (versions): length: 800/1000/1200 mm width: 410/510/610 mm height: 1600/1800/2000 mm



Dismountable Heavy load shelves



rigid construction

Pins for shelves and rails



Castors with brakes





- Stool base is made of premium stainless steel/aluminum
- Molded seats and backrests are made of integral polyurethane foam
- Color of seat and backrest covers may be chosen (accessory)



SMPS

Ø370 mm seat Seat height adjustment range: 430...630 mm Screw-type mechanism Weight max. 12 kg Ø 50 mm castors



SMP

Ø370 mm seat Seat height adjustment range: 430...630 mm Screw-type mechanism Weight max. 7,5 kg Ø 50 mm castors



SMPP

Ø370 mm seat Seat height adjustment range: 490...630 mm Hydraulic jack Weight max. 14 kg Ø 50 mm castors



SMPP-01

Ø370 mm seat
Seat height adjustment range:
490...630 mm
Hydraulic jack
Weight max. 12 kg
Ø 50 mm castors



SMPP-02

Ø370 mm seat
Seat height adjustment range:
480...670 mm
Pneumatic drive
Weight max. 13 kg
Ø 50 mm castors



SMPP-03

Ø370 mm seat Seat height adjustment range: 480...670 mm Pneumatic drive Weight max. 11 kg Ø 50 mm castors





96 | Medical furniture | 97





- Designed for convenient arrangement of instruments, medicines, and materials in medical facilities
- Made of premium stainless steel



SM-K-A

Dimensions (L/W/H): - 1060x655x915 mm Ø 125 mm castors



SM-K-A1

SI-50

Dimensions (L/W/H):

- 700x500x950 mm

Manual adjustment

900...1250 mm

SI-60-1

SI-60-2

SI-61-1

SI-61-2

Ø 80 mm castors

Dimensions (L/W/H):

1000x650x850 mm

Dimensions (L/W/H): 1200x675x850 mm Ø 80 mm castors

Dimensions (L/W/H):

1700x650x850 mm

Dimensions (L/W/H): 1800x675x850 mm Ø 80 mm castors

(without castors)

(without castors)

Height adjustment range:

Ø 125 mm castors



MEDICAL INSTRUMENT TABLES

- Designed to equip surgical departments of healthcare facilities
- Made of premium stainless steel



SI-51

Dimensions (L/W/H): - 700x500x900 mm Adjustment by foot pedal Height adjustment range: 900...1300 mm Ø 80 mm castors



SI-60

Dimensions (L/W/H): 1000x650x850 mm



Dimensions (L/W/H): 1700x650x850 mm



Dimensions (L/W/H): 1200x725x1000 mm



SI-62

Ø 150 mm castors



Dimensions (L/W/H): - 760x655x915 mm



SM-K

Push bar for convenient transportation

• Used for placing instruments, medical consumables,

MEDICAL TABLES

SM-K-01

(stainless steel) SM-K-11 (steel with polymeric coating) Dimensions (L/W/H): - 720x460x805 mm

and medicines in dressing rooms, operating rooms, and other rooms of healthcare facilities



SM-K-02

(stainless steel)

SM-K-12

(steel with polymeric coating)

Dimensions (L/W/H):

- 720x460x805 mm

Ø 80 mm castors



SM-K-03

(stainless steel)

Ø 80 mm castors

SM-K-13

(steel with polymeric coating)

Dimensions (L/W/H):

- 720x460x905 mm

Ø 80 mm castors



SM-K-04

(stainless steel)

SM-K-14

(steel with polymeric coating)

Dimensions (L/W/H):

- 720x460x905 mm

Ø 80 mm castors

SM-K-05

(stainless steel)

SM-K-15

(steel with polymeric coating)

Dimensions (L/W/H):

- 1040x700x805 mm Ø 100 mm castors



SM-K-06

(stainless steel)

SM-K-16

(steel with polymeric coating)

Dimensions (L/W/H):

- 1040x700x805 mm

Ø 100 mm castors

SM-K-07

(stainless steel)

SM-K-17

(steel with polymeric coating) Dimensions (L/W/H):

- 1040x700x905 mm

Ø 100 mm castors



SM-K-08

(stainless steel)

SM-K-18

(steel with polymeric coating)

Dimensions (L/W/H):

- 1040x700x905 mm

Ø 100 mm castors

Medical furniture | 99 98 | Medical furniture





• Made of premium stainless steel



PT-1 Medical bowl holder Single 9 L bowl, Ø500 mm base Height: 850 mm Ø 50 mm castors



PT-2
Medical bowl holder
Two 6.5 L bowls,
Ø500 mm base
Height: 850 mm
Ø 50 mm castors



PM, PIM
MEDICAL CO

MEDICAL CONTAINERS AND BUCKETS

• Made of premium stainless steel



PM-1 Waste container Dimensions (L/W/H): - 420x440x800 mm Ø 50 mm castors



PIM-1
Bucket for disposal material
Ø440 mm basis
Height: 850 mm
Ø 50 mm castors



· Made of premium stainless steel



TP-B
Trolley for linen
Dimensions (L/W/H):
- 1000x600x1000 mm
Ø 150 mm castors



Trolley for dressing boxes
Dimensions (L/W/H):
- 1200x725x1000 mm
Ø 150 mm castors

TP-SK1





- Designed for holding bottles and single-use systems with drug solutions used during therapeutic procedures in healthcare facilities
- Stand poles and metal attachments are made of premium stainless steel
- Dismountable construction with lightweight base
- Vibration and shock resistant





100 | Medical furniture | 101







SHMN-1

Height adjustment range: 0...500 mm, by means of screw locking mechanism Fixed to the wall with 2 supporting arms



SHMN-2

Height adjustment range: 0...500 mm, by means of screw locking mechanism Fixed with 1 wall guide



SHMN-3

Height adjustment range: 0...500 mm, by means of screw locking mechanism Fixed with 1 wall guide



SHMN-4

No height adjustment Fixed with 2 wall guides



Accessories



Shelf with 2 drawers



Shelf with 1 drawer

Plastic bags holder



Shelf



















Gloves and towels

dispenser holder







Clamp for wall rail

Ceiling rail and wall rail





SHMP-6

No telescopic height adjustment

With mechanical clamp and - 640x340 mm shelf



SHMP-6

No telescopic height adjustment

With mechanical clamp and - 440x340 mm shelf



SHMP-4

No telescopic height adjustment

With mechanical clamp



SHMP-3

No telescopic height adjustment

With mechanical clamp



SHMP-1

Removable, with telescopic height adjustment range above the floor: 1500...2000 mm





SHMPP-1

Suspended medical ceiling screen

Designed to separate spaces within healthcare facilities

Ceiling height: 2500...3500 mm (versions)



SH-1

Medical screen

With adjustable support and two swivel castors with brakes

The sheet is made of dense fabric with polyurethane waterproof coating Several screens can be joined together to form a unit Dimensions (L/W/H):

920x500x1760 mm

Weight: 5 kg

Medical furniture | 103 102 | Medical furniture

SA TROLLEYS FOR MEDICAL EQUIPMENT

- Used to house medical diagnostic equipment, instruments, tools, and other medical devices
- Made of engineering steel and aluminum profile coated with polymer paint
- Strong, durable construction, convenient for sanitary treatment



SA-1

Dimensions (L/W/H):

- 700x660x930 mm

Shelf dimensions (W/D):

- 550x490 mm

Height-adjustable metal shelves with working surface made of extra strong composite material

The working surface of the shelves is 6...7 mm flanged on four sides

Perforated cover made of engineering steel

Potential equalization

Accessories

Ø 125 castors

(SA-6, SA-7 Ø100 mm),

2 equipped with brakes

Handles for transportation

Shock absorbing bumpers



Endoscope holder



Screen holder and

infusion stand



Portal screen

holder

III I TII TI

11111 11111 1

.....



Lockable box





Lockable box

Lockable box

Pull-out shelf (for keyboard)

Video camera holder

SA TROLLEYS FOR MEDICAL EQUIPMENT



SA-2

Dimensions (L/W/H):

- 700x660x1200 mm

Shelf dimensions (W/D):

- 550x490 mm



SA-3

Dimensions (L/W/H):

- 700x660x1470 mm

Upper shelf dimensions (W/D): - 550x300 mm

Shelf dimensions (W/D):

550x490 mm



SA-4

Dimensions (L/W/H):

- 550x655x910 mm

Shelf dimensions (W/D):

- 490x400 mm



SA-5

Dimensions (L/W/H):

- 550x655x1180 mm

Shelf dimensions (W/D):

- 490x400 mm



SA-8

Dimensions (L/W/H):

- 1000x660x1200 mm

Shelf dimensions (W/D):

850x490 mm



SA-9

Dimensions (L/W/H):

- 1000x660x1470 mm

Upper shelf dimensions (W/D):

- 850x300 mm

Shelf dimensions (W/D):

- 850x490 mm







Dimensions (L/W/H):

- 580x620x1225 mm

Shelf dimensions (W/D):

490x400 mm

104 | Medical furniture Medical furniture | 105

TM FUNCTIONAL MEDICAL TROLLEYS



TM-1 (for procedures)

Dimensions (L/W/H):

- 790x665x1,540 mm


- Console dispenser

- Plastic bag holder

- Bottle and cups holder

The trolleys include the following:

- Designed to organize the work of medical staff
- Wide range of accessories for various purposes
- · Provide maneuverability and silent movement
- Easy to clean and disinfect





Plastic bag holder

Folding plate



Dividing inserts set



Swivel bow



Container for waste

plate



Container for collecting used sharp instruments

Towels holder



Bottle and cups holder



Basket



Infusion stand

(except TM 1, 2, 8, 9)

Endoscope holder (except TM 1, 2, 8, 9)





TM-2 (for intensive care)

Dimensions (L/W/H):

- 790x665x1.540 mm

The trolleys include the following:

- Console dispenser
- Plastic bag holder
- Bottle and cups holder
- Power filter
- Pocket for power filter cord
- Gloves and towels dispenser holder
- Container for collecting used sharp instruments
- Console of intensive care
- Key



TM-3 (for rounds)

Dimensions (L/W/H):

- 780x655x890 mm

The trolleys include the following:

- Plastic bag holder
- Basket
- Kev

TM-4 (for plastering)

Dimensions (L/W/H):

- 780x655x950 mm

The trolleys include the following:

- Plastic bag holder
- Container for collecting used sharp instruments
- Container for waste
- Towel holder
- Small working surface
- Key



TM-5 (for plastering and care)

Dimensions (L/W/H):

- 1.080x655x950 mm

The trolleys include the following:

- Plastic bag holder
- Container for collecting used sharp instruments
- Container for waste
- Towel holder
- Large working surface
- Key

TM-6 (for anesthesia)

Dimensions (L/W/H):

- 780x655x890 mm

The trolleys include the following

- Container for collecting used sharp instruments
- Container for waste
- Stainless steel plate
- Container for catheters
- Folding plate
- Key



TM-7 (for anesthesia)

Dimensions (L/W/H):

- 1,080x655x890 mm

The trolleys include the following:

- Container for collecting used sharp instruments
- Container for waste
- Stainless steel plate
- Container for catheters
- Folding plate
- Key

106 | Medical furniture Medical furniture | 107

TM FUNCTIONAL MEDICAL TROLLEYS



TM-8 (for gynecological room, instrumental)

Dimensions (L/W/H):

- 705x665x1,310 mm

The trolleys include the following:

- Swivel bowl
- Container for collecting used sharp instruments
- Waste boxes
- Container for waste
- Swivel instrumental plate
- Stainless steel plate
- Removable stainless steel cups
- Open-type bactericidal irradiator
- Key



TM-9 (for gynecological room, instrumental)

Dimensions (L/W/H):

- 705x665x1,310 mm

The trolleys include the following:

- Power filter
- Pocket for power filter cord
- Colposcope holder console



TM-10 (for emergency treatment)

Dimensions (L/W/H):

- 780x655x890 mm

The trolleys include the following:

- Container for waste
- Swivel instrumental plate
- Console of intensive care
- Infusion stand
- Oxygen bottle holder
- Power filter
- Pocket for power filter cord
- Key



TM-11 (for endoscopy)

Dimensions (L/W/H):

- 780x655x890 mm

The trolleys include the following:

- Endoscope holder stand
- Container for catheters
- Plastic bag holder
- Key





- Multipurpose couch for medical examinations and manipulations
- Durable frame
- Head section with smooth up and down adjustment by means of gas springs



Dimensions	KM-1	KM-2	KM-3	KM-4
Length, mm	1950	1950	1950	1950
Width, mm	650	650	800	800
Height, mm	650	830	650	830
Weight, kg	32	34	36	38





Upholstery colors

3096	3022	1549	1017	1089	1210	1044	6156	6279	6001	6099	5154	5141	5118	3036	7000	9011	9001





- Antistatic anti-glare glass
- Table or wall mounted
- Independent smooth brightness adjustment for each screening area
- Independent switching on of each screening area
- Low power consumption



Dimensions	NM-1	NM-2	NM-3
Length, mm	405	765	1125
Width, mm	115	115	115
Height, mm	515	515	515
Weight, kg	8	13	18

Max. brightness: (4,200...1,000) cd/m2 Single screening area size: 360x400 mm

Start-up time: 5 s

108 | Medical furniture | 109

COLOR MAP OF THE COMPANY

UPHOLSTERY COLORS*:




COLOR MAP OF THE COMPANY

UPHOLSTERY MATERIALS (KMM, KMMP)*:



COATING COLORS FOR METALLIC SURFACES*:



^{*} ATTENTION! Real-life colors may differ: graphic prints or screen images do not provide exact match.

110 | Color map | 111

^{*} ATTENTION! Real-life colors may differ: graphic prints or screen images do not provide exact match.

www.medin.by

For more information, please visit our website