



SHDV, SHDVP

MEDICAL STANDS FOR LONG-LASTING INFUSIONS



LEARN MORE

- Designed for holding bottles and single-use systems with drug solutions used during therapeutic procedures in healthcare facilities
- Stand poles and metal attachments are made of premium stainless steel
- Dismountable construction with lightweight base
- Vibration and shock resistant



SHDV-1

Height: 1855 mm
No height adjustment
No castors



SHDVP-1

Height: 1900 mm
No height adjustment
Ø 50 mm castors



SHDVP-3

Height adjustment range:
1,500...2325 mm, by means
of automatic locking
Z-base
Ø 80 mm castors



SHDV-2

Height adjustment range:
1400...2200 mm, by means
of screw locking mechanism
No castors



SHDVP-2

Height adjustment range:
1430...2230 mm, by means
of screw locking mechanism
Ø 50 mm castors



SHDVP-4

Height adjustment range:
1525...2340 mm, by means
of automatic locking
Ø 80 mm castors

SHMN MEDICAL WALL STANDS



SHMN-1
Height adjustment range: 0...500 mm,
by means of screw locking mechanism
Fixed to the wall with 2 supporting arms



SHMN-2
Height adjustment range: 0...500 mm,
by means of screw locking mechanism
Fixed with 1 wall guide



SHMN-3
Height adjustment range: 0...500 mm,
by means of screw locking mechanism
Fixed with 1 wall guide



SHMN-4
No height adjustment
Fixed with 2 wall guides

Accessories



Shelf with 2 drawers



Shelf with 1 drawer



Shelf



Ceiling rail and wall rail



Swivel bowl



Plastic bags holder



Gloves and towels
dispenser holder



Bottle holder



Clamp for wall rail



Equipment holder

SHMP MEDICAL CEILING STANDS



SHMP-6
No telescopic height
adjustment
With mechanical clamp and
- 640x340 mm shelf



SHMP-6
No telescopic height
adjustment
With mechanical clamp and
- 440x340 mm shelf



SHMP-4
No telescopic height
adjustment
With mechanical clamp



SHMP-3
No telescopic height
adjustment
With mechanical clamp



SHMP-1
Removable, with
telescopic height
adjustment range
above the floor:
1500...2000 mm